Kvalitetskontroll av resistensbestemmelse

AFA- kurs januar 2026

Alexander Husby

Avdeling for medisinsk mikrobiologi, St. Olavs hospital

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The European Committee on Antimicrobial Susceptibility Testing

NordicAST Nordic Committee on Antimicrobial Susceptibility Testing

EUCAST is a scientific committee of European breakpoint committees. EUCAST provides guidance documents to aid clinical microbiology laboratories in the interpretation of antimicrobial susceptibility test results.

Further details about EUCAST and its work can be found on the EUCAST website.

Most accessed links

NordicAST	Organisation	Methods	Breakpoints	Quality control	Resistenspaneler	FAQ	Workshops	News
Links	RSS							

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
- New "Routine and extended internal quality control..." 19 hours ago
- New publication on how to interpret antifungal age... 22 hours ago
- Version 16.0 of the breakpoint tables: what's new 4 days ago
- Minor change to the EUCAST Screening method for th... 20 days ago
- Webinar to update on version 16.0 Breakpoint Table... 20 days ago
- Consultation on minor changes to version 16.0 Brea...

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MIC determination (broth microdilution according to ISO standard 20776-1 except for mecillinam and fosfomycin where agar dilution is used)

Medium: Mueller-Hinton broth

Inoculum: 5×10^5 CFU/mL

Incubation: Sealed panels, air, $35 \pm 1^\circ\text{C}$, $18 \pm 2\text{h}$

Reading: Unless otherwise stated, read MICs at the lowest concentration of the agent that completely inhibits visible growth.

Quality control: *Escherichia coli* ATCC 25922. For agents not covered by this strain and for control of the inhibitor component of beta-lactam inhibitor combinations, see EUCAST QC Tables.

MIC determination (broth microdilution according to ISO standard 20776-1)

Medium: Cation-adjusted Mueller-Hinton broth + 5% lysed horse blood and 20 mg/L β -NAD (MH-F broth)

Inoculum: 5×10^5 CFU/mL

Incubation: Sealed panels, air, $35 \pm 1^\circ\text{C}$, $18 \pm 2\text{h}$

Reading: Unless otherwise stated, read MICs at the lowest concentration of the agent that completely inhibits visible growth. See "EUCAST Reading Guide for broth microdilution" for further information.

Quality control: *Haemophilus influenzae* ATCC 49766. For agents not covered by this strain and for control of the inhibitor component of beta-lactam inhibitor combinations, see EUCAST QC Tables.

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



- Egne BP tabeller for RAST ligger på EUCAST.

Escherichia coli

Zone diameter breakpoints for RAST directly from blood culture bottles

EUCAST RAST Breakpoint Tables v. 6.1, valid from 2023-06-07

EUCAST rapid disk diffusion method directly from positive blood culture bottles
 Medium: Mueller-Hinton (MH) agar
 Inoculum: 125±25 µL directly from a positive blood culture bottle
 Incubation: Air, 35±1°C
 Incubation time: 4, 6, 8 and 16-20 hours
 General reading instructions: Inhibition zones should only be read when the growth is confluent and zone edges are clearly visible.
 Reading 4, 6 and 8 hours: Remove the lid and read zone diameters from the front against a dark background illuminated with reflected light.
 Reading 16-20 hours: Read zone diameters from the back of the plate against a dark background illuminated with reflected light.
[QC for implementation of RAST](#)

Antimicrobial agent	Disk content (µg)	4 hours			6 hours			8 hours			16-20 hours		
		S ≥	ATU	R <	S ≥	ATU	R <	S ≥	ATU	R <	S ≥	ATU	R <
Ampicillin	10	12	10-11	10	13	11-12	11	13	11-12	11	14	12-13	12
Amoxicillin-clavulanic acid	20-10	16	14-15	14	17	15-16	15	17	15-16	15	19	17-18	17
Piperacillin-tazobactam	30-6	17	14-16	14	18	15-17	15	18	15-17	15	17	15-16	15
Temocillin	30	50	13-14	13	50	15-16	15	50	15-16	15	50	17-18	17
Cefotaxime ¹	5	15	13-14	13	16	14-15	14	17	15-16	15	16	14-15	14
Ceftazidime ¹	10	15	12-14	12	16	14-15	14	17	15-16	15	17	15-16	15
Ceftazidime-avibactam	10-4	12	10-11	10	12	10-11	10	12	10-11	10	13	11-12	11
Ceftolozane-tazobactam	30-10	16	14-15	14	18	16-17	16	18	16-17	16	20	16-19	16
Imipenem ²	10	16	12-15	12	17	13-16	13	17	13-16	13	17	12-16	12
Imipenem-relebactam	10-25	13	11-12	11	15	13-14	13	15	13-14	13	17	15-16	15
Meropenem ²	10	17	15-16	15	17	15-16	15	17	15-16	15	15	13-14	13
Meropenem-vaborbactam	20-10	16	12-15	12	16	14-15	14	17	15-16	15	16	14-15	14
Ciprofloxacin	5	17	14-16	14	19	16-18	16	20	17-19	17	21	18-20	18
Levofloxacin	5	16	14-15	14	18	15-17	15	18	15-17	15	23	15-22	15
Amikacin ³	30	(15)	(13-14)	(13)	(15)	(13-14)	(13)	(15)	(13-14)	(13)	(13)	(11-12)	(11)
Gentamicin ³	10	(14)	(12-13)	(12)	(14)	(12-13)	(12)	(14)	(12-13)	(12)	(13)	(11-12)	(11)
Tobramycin ³	10	(14)	(12-13)	(12)	(15)	(13-14)	(13)	(15)	(13-14)	(13)	(13)	(11-12)	(11)
Trimethoprim-sulfamethoxazole	1.25-23.75	12	10-11	10	14	12-13	12	14	12-13	12	14	12-13	12

Notes

1. Cephalosporin breakpoints for *E. coli* will detect all clinically important resistance mechanisms. The presence or absence of an ESBL does not in itself influence the categorisation of susceptibility. However, ESBL detection and characterisation are recommended for public health and infection control purposes.

See document [EUCAST RAST screening for resistance mechanisms \(link to document\) for screening cut-offs](#).

2. Carbapenem breakpoints for *E. coli* will detect all clinically important resistance mechanisms. The presence or absence of a carbapenemase does not in itself influence the categorisation of susceptibility. However, carbapenemase detection and characterisation are recommended for public health and infection control purposes.

See document [EUCAST RAST screening for resistance mechanisms \(link to document\) for screening cut-offs](#).


3. Aminoglycoside breakpoints distinguish between isolates without and with resistance mechanisms. For blood stream infections, EUCAST recommends that aminoglycosides are used in combination with other active therapy.

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European Committee on Antimicrobial Susceptibility Testing

Routine and extended internal quality control
for MIC determination and disk diffusion
as recommended by EUCAST

Version 14.0, valid from 2024-01-01

This document should be cited as

"The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. Version 14.0, 2024. <http://www.eucast.org>."

General	Page
Notes	1
Changes	2

Routine quality control	Page
Recommended strains for routine quality control	4
<i>Escherichia coli</i> ATCC 25922	6
<i>Pseudomonas aeruginosa</i> ATCC 27853	8
<i>Staphylococcus aureus</i> ATCC 29213	9
<i>Enterococcus faecalis</i> ATCC 29212	11
<i>Streptococcus pneumoniae</i> ATCC 49619	12
<i>Haemophilus influenzae</i> ATCC 49766	14
<i>Campylobacter jejuni</i> ATCC 33560	15
<i>Mannheimia haemolytica</i> ATCC 33396	16
<i>Bacteroides fragilis</i> ATCC 25285	17
<i>Clostridium perfringens</i> ATCC 13124	18
Control of the anaerobic environment	19
Control of the inhibitor component of β -lactam-inhibitor combinations	20

Extended quality control for detection of resistance mechanisms with disk diffusion	Page
ESBL production in <i>Enterobacterales</i>	23
Methicillin resistance in <i>Staphylococcus aureus</i>	23
<i>vanB</i> -mediated glycopeptide resistance in enterococci	23
Acquired aminoglycoside-modifying enzyme (high-level aminoglycoside resistance) in enterococci	24
Reduced susceptibility to β -lactam agents due to PBP mutations in <i>Haemophilus influenzae</i>	25

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
UK NEQAS Microbiology		Antimicrobial susceptibility	Laboratory : 1487
		Distribution : 4297	Page 1 of 11
		Dispatch Date : 04-Jun-2018	
Intended Result	Your Report		Your Score
Specimen 4430	Streptococcus pneumoniae from cerebrospinal fluid		
Cefotaxime	intermediate	intermediate	Not scored
Meningitis/cefotaxime	intermediate	not examined	Not scored
Ceftriaxone	intermediate	not examined	Not scored
Meningitis/ceftriaxone	intermediate	not examined	Not scored
Penicillin	resistant	intermediate	Not scored
Meningitis/PEN	resistant	resistant	2
Vancomycin	susceptible	not examined	Not scored
Specimen 4431	Klebsiella pneumoniae from blood		
Amikacin	susceptible	not examined	Not scored
Amoxicillin	resistant	not examined	Not scored
Amoxicillin-clavulanic acid	resistant	resistant	2
Ampicillin	resistant	resistant	2
Cefotaxime	resistant	resistant	2
Cefazidime	resistant	resistant	2
Ceftriaxone	resistant	not examined	Not scored
Cefuroxime	resistant	resistant	2
Ciprofloxacin	susceptible	susceptible	2
Colistin	susceptible	not examined	Not scored
Ertapenem	resistant	resistant	2
Gentamicin	resistant	resistant	2
Imipenem	intermediate	susceptible	0
Meropenem	resistant	intermediate	1
Piperacillin-tazobactam	resistant	resistant	2
Tigecycline	susceptible	susceptible	2
Tobramycin	resistant	resistant	2
ESBL	negative	not examined	Not scored
ArgC	negative	not examined	Not scored
Carbapenemase	positive	positive	2

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